




6505 013723425 SODIUM CHLORIDE INJ PG I.V Bag 500ML 2 each

(lactated ringers may be substituted for SODIUM CHLORIDE)



6510 000583047 BANDAGE GAU4-1/2"100S PG Gauze Kerlix



6510 001594883 DRESSING FIRST AID EAArmy First-aid



0 002011755 BANDAGE 37X37X52IN EA

Cravats, green

(Green or khaki is acceptable)



6510 002017425 DRESS FLD 11-3/4 IN EA Abd dressing



6510 007863736 PAD ISOPROPYL ALCOHOL PG square alcohol pads



6510 009137909 BAND ADH .75X3 IN300S BX

Band-Aids



6510 009268882 ADHESIVE TAPE SURG 1" PG


white

one r



6510 010100307 PAD POV-IOD IMPRE100S PG

square pov-iod pads



6510 014640826 SPONGE SURG 2X2 100S PG

white-gauze



6515 002267692 GLOVE EXAM LARGE 100S PG non-sterile




6515 009582232 AIRWAY PHARYNX 80MM12S BX J-tube(s-adult) 1 each

PHARYNX

100MM BX J-tube(adult)


(The size is indicated near the flat tip)
each



6515 009357138 SCISSORS BANDAGE 7.25" EA

1 each

(Either type of scissor is acceptable)



6515 011885316 TUBE DRAIN 1X18" 200S PG constricting band 6 each
(there are several types of constricting bands available; as long as there is at least one)

SPLINT, UNIVERSAL
6515-01-225-4681
BRENNER METAL PRODUCTS

UNROLL THE SPLINT AND IT CAN BE USED
AS FOLLOWS:

1. THE SPLINT MAY BE USED IN A
SIMPLE "V" SHAPE.



2. BY BENDING IT TO A "U" SHAPE,
YOU'RE STRENGTHENING
SPLINT.



3. THE SPLINT CAN THEN BE
APPLIED TO THE WRIST.
FOREARM, UPPER ARM, LOWER





6515 013156227 CATHETER/NDL 18GA 50S PG I.V Needle

2 each

absolutely nothing larger than 18 gauge and two inches long; the smaller the number the bigger the needle.



6545 009129870 CASE #3 10X4-1/2X18IN EA

Combat Lifesaver Bag

(an a5.56 ammo can is an excellent substitute)